

Mandatory Medical Device Problem Reporting Form for Industry

CANADA VIGILANCE - MEDICAL DEVICE PROBLEM REPORTING PROGRAM

If more space is required, please attach additional sheets

Fields required to be completed for updates/final reports are indicated by an $^{\star}\,$

☐ Preliminary ☐ Update ☒ Final ☐ Preliminary & Final If "preliminary" only, anticipated date for the final report:	
Y-MM-DD)	
Health	
Canada:	
Y-MM-DD)	
Y-MM-DD)	
A	

Canada Canada

E. INVESTIGATION INFORMATION

1. Trade/Brand Name *	1. Investigative Actions and Timeline
SYNCHROMED II	Text for this section follows on the subsequent page(s)
2. Control/Lot/Serial No.	
NGV477565H	
3. Expiration Date	
2014-10-28	
4. i. Device Classification	
ii. Device License No.	
63074	
iii. Device Identification No	
426043 iv. Manufacturer's Medical Device Identifier (catalogue/model no.)	
8637-40	
5. Software Version	
Not Applicable	
6. Age of Device	
Manufacture Date:2013-05-02	This section only applies for preliminary & final, and final reports
7. How long was the device in use?	2. Root Cause of Problem
Implant Date:2014-02-24	The root cause of the motor stall was determined to be related to residue on the
8. Was the device labelled as sterile?	gear shaft of the motor gear train.
⊠ Yes □ No	
9. Availability of device for evaluation	
☐ Destroyed ☐ Returned to Manufacturer/Importer	
☐ Neither (with explanation)	
D. COMPLAINANT INFORMATION	
1. Complainant is a:	
☐ Consumer ☐ Health professional ☐ Other	
2. Name of Complainant	
Femir Nesib	
3. Name of Health Care Facility (if applicable)	
Health Sciences Centre Winnipeg	
	3. Corrective Actions taken as a result of the investigation
4. Address	
820 Sherbrook St, Winnipeg, MB R3A 1R9, Canada	No corrective actions have been taken as a result of this event. This failure mode will continue to be monitored and trended.
	will continue to be morniored and deflued.
5. Telephone No. and/or E-mail Address	
12047872365	
Privacy Notice Statement: For the purposes of the Canada Vigilance - Medical Device Problem Reporting Program, information related to the identity	
of the complainant and/or reporter will be protected as personal information	
under the <i>Privacy Act</i> , under the <i>Access to Information Act</i> in the case of an access to information request. For details with regard to personal	
information collected under this program, visit the Personal Information Bank;	
Health Canada; Health Products and Food Branch; Branch Incident Reporting	

C. MEDICAL DEVICE INFORMATION

System; HC PPU 088 at: http://infosource.gc.ca/inst/1476/1476-fedemp00-

eng.asp

Disclaimer

Medtronic, Inc. (Medtronic) is submitting this report to comply with regulatory reporting requirements. This report is based upon information obtained by Medtronic, which the company may not have been able to fully investigate or verify prior to the date the report was required to be submitted. Medtronic has made reasonable efforts to obtain more complete information in the time allotted and has provided as much information as is available to the company as of the submission date of this report. This report does not constitute an admission or a conclusion by the regulatory authority, Medtronic, or its employees that the device, Medtronic or its employees caused or contributed to the event described in the report. In particular, this report does not constitute an admission by anyone that the product described in this report has any "defects" or has "malfunctioned". Medtronic objects to the use of thesewords and others like it because of the lack of definition and the connotations implied by these terms. This statement should be included with any information or report disclosed to the public under applicable disclosure laws.

B. 5. Details of Incident

It was reported that the patient heard their alarm on 2020-APR-07 and went into the clinic to investigate the alarm. Upon interrogation, it read motor stall occurred on 2020-APR-07 at 8:40 PM and recovered on 2020-APR-08 at 11:59 AM. The patient did not have an MRI, and it was unsure if the patient had any EMI exposure. There were no reported patient symptoms. Additional information was received in response to a request for follow-up. It was reported interrogating the pump a few times helped to resolve the issue and it was confirmed by the pump message "motor stall recovery" had occurred. It was confirmed that the patient denied going for any diagnostic testing such as an MRI, but was noted that the patient worked with computers and digital technologies. The HCP had given the patient a prescription for oral baclofen in case the issue occurred again and the patient was encouraged to call the clinic with any concerns. Additional information was received. No replacement was planned at this time.

Additional information was received on 2020-APR-27. It was indicated the pump stalled again and had not recovered. There was no EMI or magnetic interaction (MRI).

Further information received indicated the pump was replaced on 2020-MAY-05 and the patient had recovered and was doing well.

Further information received on 2020-AUG-14 indicated that logs showed tube set occurred.

E. 1. Investigative Actions and Timeline

It was reported that multiple motor stalls occurred, resulting in tube set. As the event was reportable to a regulatory authority and indicated a possible failure of a device, labeling, or packaging to meet any of its specifications, an investigation was required. Additional information was required to investigate the event and/or determine the cause of the event. The healthcare provider (HCP) was contacted via a manufacturer representative to investigate contributing factors. The HCP's response via the manufacturer representative indicated that the patient denied going for any diagnostic testing i.e. MRI. However, the patient said he works with computers, and digital technologies. The device was returned as part of this investigation. Analysis found residue in the motor gear train. Analysis identified wearing on the upper shaft of gear number two. There was no indication that the event was related to a possible manufacturing issue, so a Device History Record review was not performed. The investigation determined that this event was similar to an event that had already been investigated, and another investigation is not necessary. See PUMP MOTOR CORROSION AND WEAR. Analysis found residue in the motor gear train. Analysis identified wearing on the upper shaft of gear number two. The root cause, if available, is documented in the referenced investigation record.

Salazar, Glenda (HC/SC)

From: Pascal, Laura < laura.d.pascal@medtronic.com>

Sent: 2020-09-23 10:47 AM **To:** mdpr / dimm (HC/SC)

Subject: MDPR - FINAL - Medtronic file # 703703023 - Date of submission: 23-Sep-2020

Attachments: 703703023 - MDR.pdf

Please find attached a FINAL – 30 DAY – MDPR. Medtronic File no. 703703023

Best regards,

Laura Pascal

Compliance/Audit Specialist CANADA

Medtronic

16720 Transcanadienne, Suite 200 Kirkland, Québec H9H 4M7 Bureau/Office: 514-693-7146 | Fax: 514-693-7278

buleau/Office. 514-095-7140 | 1 ax. 514-095-7278

laura.d.pascal@medtronic.com



[CONFIDENTIALITY AND PRIVACY NOTICE] Information transmitted by this email is proprietary to Medtronic and is intended for use only by the individual or entity to which it is addressed, and may contain information that is private, privileged, confidential or exempt from disclosure under applicable law. If you are not the intended recipient or it appears that this mail has been forwarded to you without proper authority, you are notified that any use or dissemination of this information in any manner is strictly prohibited. In such cases, please delete this mail from your records. To view this notice in other languages you can either select the following link or manually copy and paste the link into the address bar of a web browser: http://emaildisclaimer.medtronic.com