

Mandatory Medical Device Problem Reporting Form for Industry

CANADA VIGILANCE - MEDICAL DEVICE PROBLEM REPORTING PROGRAM

If more space is required, please attach additional sheets

Fields required to be completed for updates/final reports are indicated by an $^{\star}\,$

			Page <u>1</u> of <u>3</u>
A. REPORTER INFORMATI	ON		
1. i. Reporter Type		3. Reporter File No. *	
☐ Manufacturer ☐ Importer		0703667127	
In the case where the reporter is the importer:		4. Health Canada File No. (if applicable) *	
ii. Did the importer report the incident to the manufacturer?			
∑ Yes □ No		5. Type of Report *	
iii. Is the importer also submitting the report on behalf of the manufacturer?		☐ Preliminary ☑ Update ☐ Final ☐ Preliminary & Final If "preliminary" only, anticipated date for the final report:	
Yes	☐ No		(YYYY-MM-DD)
2. Reporter Contact Information *		If "update/final", date the previous report was submitted to Health Canada:	
Laura Pascal - RA/QA Dept. Medtronic Canada Phone: 514-693-7146			
Laura.d.pascal@medtronic.com		2020-10-27	(YYYY-MM-DD)
	Manufacturer		Importer
7. Name and Address	MEDTRONIC NAVIGATION, INC. 826 COAL CREEK CIRCLE, LOUISVILLE, CO, US Postcode:80027 Tel:(+1-720)8903200 Fax:(+1-720)8903500		MEDTRONIC CANADA ULC 99 HEREFORD STREET, BRAMPTON, ON, CANADA Postcode:L6Y 0R3 Tel:(800)268-5346 Fax:(905) 460-3992
8. Health Canada assigned company identification number (if	114892		106916
9. Establishment License Number (if applicable):	NA		35
B. INCIDENT INFORMATIO	N		
1. Classification of Incident *		5. Details of Incident	
i. □ 10-Day ⊠ 30-Day ii. ☑ Canadian □ Foreign iii. □ Investigational testing □ Special Access Program □ Radiation emitting device (if applicable)		Medtronic received information regarding a navigation system being used in a cranial biopsy. It was reported that the site couldn't get the system to recognize the emitter box. The system could still be used Clinically as optical was working, but was tagged as such until the box could be tested/replaced. No patient was present at the time of the event.	
2. Date of Incident		Concomitant Medical Products:	
2020-03-13		Product ID: 9733597, Product Type: CABLE 9733597 EXTERNAL AXIEM PWR/DATA.	
3. Reporter's Awareness Date			
2020-03-13			
4. Patient Consequences			
No patient consequences were reported			



E. INVESTIGATION INFORMATION

1. Trade/Brand Name *	1. Investigative Actions and Timeline
STEALTHSTATION® S7'	Text for this section follows on the subsequent page(s)
2. Control/Lot/Serial No.	
4501123545	
3. Expiration Date	
4. i. Device Classification	
□ I □ II ⊠ III □ IV ii. Device License No.	
79099	
iii. Device Identification No	
iv. Manufacturer's Medical Device Identifier (catalogue/model no.)	
9733857 5. Software Version	
Not Applicable 6. Age of Device	
7. How long was the device in use?	This section only applies for preliminary & final, and final reports 2. Root Cause of Problem
UNKNOWN 8. Was the device labelled as sterile?	The lemo connector on the returned cable had been previously opened and the alignment sleeve is missing. Without the sleeve, the pin block is not set will not
☐ Yes ☐ No	connect correctly. The cable is not usable as is.
9. Availability of device for evaluation	1
☐ Destroyed ☐ Returned to Manufacturer/Importer	
Neither (with explanation)	
Still in Use	
Product ID: 9733597, Product Type: CABLE 9733597 EXTERNAL AXIEM	
D. COMPLAINANT INFORMATION	
1. Complainant is a:	
☐ Consumer ☐ Health professional ☐ Other	
2. Name of Complainant	
Unknown Healthcare Professional	
3. Name of Health Care Facility (if applicable)	
Foothills Medical Centre	
	3. Corrective Actions taken as a result of the investigation
4. Address	No corrective actions have been taken as a result of this event. We will continue
1403-29th Street NW Calgary AB T2N 2T9	to monitor for trends as more definitive data is collected.
5. Telephone No. and/or E-mail Address	
Tel#: +1UNKNOWN#	
Privacy Notice Statement: For the purposes of the Canada Vigilance - Medical Device Problem Reporting Program, information related to the identity of the complainant and/or reporter will be protected as personal information under the <i>Privacy Act</i> , under the <i>Access to Information Act</i> in the case of an access to information request. For details with regard to personal information collected under this program, visit the Personal Information Bank;	

C. MEDICAL DEVICE INFORMATION

System;HC PPU 088 at: http://infosource.gc.ca/inst/1476/1476-fedemp00-

eng.asp

Disclaimer

Medtronic, Inc. (Medtronic) is submitting this report to comply with regulatory reporting requirements. This report is based upon information obtained by Medtronic, which the company may not have been able to fully investigate or verify prior to the date the report was required to be submitted. Medtronic has made reasonable efforts to obtain more complete information in the time allotted and has provided as much information as is available to the company as of the submission date of this report. This report does not constitute an admission or a conclusion by the regulatory authority, Medtronic, or its employees that the device, Medtronic or its employees caused or contributed to the event described in the report. In particular, this report does not constitute an admission by anyone that the product described in this report has any "defects" or has "malfunctioned". Medtronic objects to the use of thesewords and others like it because of the lack of definition and the connotations implied by these terms. This statement should be included with any information or report disclosed to the public under applicable disclosure laws.

E. 1. Investigative Actions and Timeline

The CABLE 9733597 has been received at the USA Global Business Unit and investigation was re-opened.

Investigation Summary - It was reported that the site couldn't get the system to recognize the emitter box. The system could still be used Clinically as optical was working but was tagged as such until the box could be tested/replaced. No patient was present at the time of the event. As the event was reportable to a regulatory authority and indicated a possible failure of a device, labeling, or packaging to meet any of its specifications, an investigation was required. The event was investigated to determine cause. The system was serviced in the field and the external axiem cable was replaced.

Hardware analysis determined that the lemo connector on the returned cable had been previously opened and the alignment sleeve is missing. Without the sleeve, the pin block is not set will not connect correctly. The cable is not usable as is. As the reported issue was unrelated to a software issue, a software analysis was not performed. There was no indication that the event was related to a possible manufacturing issue, so a Device History Record review was not performed. The investigation determined that this event was similar to an event that had already been investigated, and another investigation is not necessary. See Test Track Number 2257. the in-house analysis indicates the event is similar to events investigated per Test Track # 2257.

Aijaz, Naeem (HC/SC)

From: Pascal, Laura < laura.d.pascal@medtronic.com>

Sent: 2020-10-27 11:42 AM **To:** mdpr / dimm (HC/SC)

Subject: MDPR - UPDATE - Medtronic file #703667127 - Date of submission: 27-Oct-2020

Attachments: 703667127 - MDR.pdf

Please find attached an Update - 30 DAY - MDPR. Medtronic File no. 703667127

Best regards,

Laura Pascal

Compliance/Audit Specialist CANADA

Medtronic

16720 Transcanadienne, Suite 200 Kirkland, Québec H9H 4M7 Bureau/Office: 514-693-7146 | Fax: 514-693-7278

laura.d.pascal@medtronic.com



[CONFIDENTIALITY AND PRIVACY NOTICE] Information transmitted by this email is proprietary to Medtronic and is intended for use only by the individual or entity to which it is addressed, and may contain information that is private, privileged, confidential or exempt from disclosure under applicable law. If you are not the intended recipient or it appears that this mail has been forwarded to you without proper authority, you are notified that any use or dissemination of this information in any manner is strictly prohibited. In such cases, please delete this mail from your records. To view this notice in other languages you can either select the following link or manually copy and paste the link into the address bar of a web browser: http://emaildisclaimer.medtronic.com